

Aflatoxin M1 ULTRA ELISA Quantitative

For the quantitative detection of Aflatoxin M1 in milk, skim milk powder, and yogurt.

Kit Components

- 96 Antibody-coated Microwells
- Mixing wells
- Aflatoxin M1 Standards
- Ready to Use Conjugate
- M1 Free Skim Milk
- Substrate Reagent
- Stop Solution
- PBS-T Powder




Catalog No: 961AFLM01C-ULTRA

Required Equipment Not Supplied with Kit


- Single or multi-channel pipettor with 10, 100, 200 and 1000 μ L tips
- Microtubes
- Wash bottle
- Absorbent paper towels
- Centrifuge (and tubes)
- Microplate reader with 450 nm filter
- Triton X-100 (20% in deionized water, v/v)
- Yogurt diluent for yogurt sample (Cat# 937YOG001)

Reagents Provided


1x Pouch	Antibody coated microwell plate		96 wells (12 eight-well strips) in a microwell holder coated with a mouse anti-aflatoxin monoclonal antibody, <i>Ready-To-Use</i>
1x Plate	Mixing wells	Red	96 non-coated wells (12 eight-well strips) in a microwell holder, <i>Ready-To-Use</i>
6x Vials	Aflatoxin M1 standards	Black cap	8.0 mL/vial of Aflatoxin M1 at the following concentrations: 0.0, 5.0, 15.0, 50.0, 150.0, 500.0 pg/mL (ppt), <i>Ready-To-Use</i>
1x Bottle	Aflatoxin HRP-conjugate	Green cap	12 mL of Aflatoxin conjugated to horseradish peroxidase in buffer with preservative, <i>Ready-To-Use</i>
1x Bottle	Substrate reagent	Blue cap	12 mL stabilized tetramethylbenzidine (TMB), <i>Ready-To-Use</i>
1x Bottle	Stop solution	Red cap	12 mL Acidic Solution, <i>Ready-To-Use</i>
1x Pouch	Washing buffer		PBS with 0.05% Tween20®, bring to 1 liter with distilled water and store refrigerated.
1x Bottle	M1 free skim milk	White cap	12 mL of skim milk, <i>Ready-To-Use</i>



1 Transfer 1.2 mL of each standard and sample into microtubes. If running singlets, scale the volume down accordingly.




2 Use a multichannel pipettor to transfer 200 μ L aliquots of standards and samples from the microtubes into the Antibody Coated Wells and seal plate. Incubate for 20 minutes as the first incubation.



3 Empty the contents from the Antibody Coated Wells and discard into basin. Wash the wells by filling with the reconstituted PBS-Tween wash buffer, empty into basin, and repeat for a total of 3 washings. Tap the wells (face down) on a layer of absorbent paper.

4 Repeat step 2. Transfer 200 μ L aliquots of standards and samples. Incubate for 20 minutes as the second incubation. Refer to photos and instructions in step 2.




5 During the second incubation, dispense 150 μ L of standard or sample from the microtube into the red mixing wells provided by the kit and add 150 μ L of the conjugate to each red mixing well. Mix by priming pipettor at least 3 times. If running singlets, scale the volume down accordingly. *Note: Operator must record the location of each Standard and Sample throughout test.*

6 After the second incubation, wash the plate by repeating step 3. Refer to photos and instructions in step 3.




7 Transfer 100 μ L of the conjugate mixture from each mixing well in step 5 to a corresponding Antibody Coated well. Seal and incubate for 20 minutes.

8 Repeat step 3. Refer to photos and instructions in step 3 for a total of **5 washings**.



9 Add 100 μ L of enzyme substrate (TMB) to each Antibody Coated well and incubate for 10 minutes. Cover to avoid direct light. A blue color will develop.

10 Use a multi-channel pipettor to transfer 100 μ L of stop solution to the Antibody Coated wells. The blue color will change to yellow.



11 Read the optical density of each microwell with a microplate reader at 450 nm using an air blank or a differential filter of 630 nm.