

Cost effective cleaning validation for allergens

Effective cleaning is usually identified as a pre-requisite for most GMP and HACCP plans in the food industry and cleaning is often considered a critical control point (CCP) for allergen control. Cleaning is designed to remove food residues that contain many common components such as ATP (adenosine triphosphate), protein and sugars. Some of these foods may also contain allergens. The more effective the cleaning procedure, then the lower the amount of food residues, and hence the lower the risk. Using the most sensitive detection methods gives the greatest assurance of cleaning efficacy.

The principle of broad spectrum monitoring methods together with specific detection methods is well established in monitoring and managing risk. In microbiology analysis for example, total bacteria counts, coliforms and *Listeria* spp are used as overall indicators and then specific tests for specific pathogens are used as required. Here we describe a similar combined approach of pre-validation and monitoring of cleaning for allergen management using a combination of three highly sensitive detection methods.

Detection methods

ATP bioluminescence provides an immediate direct objective test of cleaning efficacy well established for >30 years and detects a very broad range of foodstuffs. Recent developments in ATP bioluminescence have improved detection capabilities and sensitivity and at this level it is capable of detecting food residues below the limit of detection of specific allergen tests.

The new EnSURE instrument and SuperSnap reagent swab (Hygiena) provide additional sensitivity with low background noise and low variation for precise, accurate and consistent results. This means the EnSURE is x10 more sensitive than Hygiena SystemSURE Plus with UltraSnap swabs and x100 more sensitive than other ATP systems (see Table 1; for further details contact



Hygiena for a copy of the Comparative report). The results are quantitative and give a linear response to increasing amounts of food residue. SuperSnap also provides more robustness and tolerance to harsh materials at extremes of pH and in the presence of sanitiser, e.g. it is not affected by 1000ppm hypochlorite.

Most allergens are glycoproteins and can be detected by protein tests (such as the biuret method, ALLER-Snap) however, this non-specific protein test cannot differentiate non-allergenic protein

from true allergens. Protein tests can detect allergenic foodstuffs, and for maximum sensitivity (1–3µg protein) the test is run at elevated time and temperature combinations such as 37°C for 30min. The results are semi-quantitative and the scope and sensitivity of the protein test is 10–100 ppm for certain allergenic foods.

Specific allergen tests such as lateral flow formats were originally designed to detect the presence of the allergen in foodstuffs and certain extraction procedures are

The New Rapid Test Platform For The 21st Century

EnSURE

The **EnSURE** system measures a novel bioluminogenic test reaction that generates light when enzymes that are characteristic of specific bacteria react with specialised substrates to produce light. The **EnSURE** system quantifies the light into an objective value that is displayed and stored in its memory.

3 NOVEL TESTS

Micro-Snap


A rapid test for enumeration of Enterobacteriaceae and E.coli. Detects single figure organisms in **JUST 7 HOURS**

Zymo-Snap

A 3 minute test of industrial important enzymes such as Alkaline Phosphatase and Proteases

SuperSnap

Highly sensitive hygiene monitoring for more demanding cleaning verification such as allergen control



hygiena
INTERNATIONAL

Call Now On:
+44 (0)1923 818821
www.hygiena.net

required for optimal performance. This technology has been extended for surface hygiene testing for cleaning validation where the limit of detection is claimed to be 1–20µg or ppm, however, studies have shown that they only achieve 4–27% recovery and in practice give a qualitative presence/absence result.

Table 2 demonstrates results from a factory trial where high sensitivity ATP and high sensitivity protein tests provided an effective monitoring tool as part of the allergen management programme. Before cleaning all test results were positive and after cleaning most were negative. The ATP test detects residues below that of protein tests and specific allergens were not detected, thus confirming that the highest level of cleaning had been achieved and allergens were absent.

Detecting allergens in a ready meal factory

A production facility manufacturing ready meals and vegetable dishes for major supermarket retailers also makes a nut product on a less frequent basis. The site needed to be sure that its cleaning had been effective to remove nut allergens after the manufacturing of nut

Table 1: Comparative sensitivity of new ATP systems

Parameter	ATP test systems		
	Hygiena Pi 102 and SuperSnap	EnSURE and SuperSnap	Others
Sensitivity (Limit of Detection) (fmol ATP)	0.01	0.1	1.0 to 10.0
Repeatability (CV%)	12%	9%	26 to 123%

products and before releasing the production area back to general manufacturing. The products contained three different tree nuts but, for the sake of completeness nine nut allergens were tested in the cleaning validation exercise and all nine nut allergens needed to be shown to be absent before release of the lines and equipment.

An off-site contract laboratory was used to conduct specific Elisa-based allergen tests with a turnaround time of 10 working days during which the production facility could not be used, thus losing valuable production time. A minimum of ten different samples were taken at various points of the production facility and each sample was tested for nine tree nut allergens at considerable cost.

Previous cleaning validation

exercises using only the specific allergen tests had not always passed first time, thus requiring repeat testing and the production line out of use for further 10 days. This was an extremely costly exercise, and the facility needed a faster, more reliable and cost effective way to validate the cleaning.

The EnSURE luminometer with SuperSnap gives a high sensitivity ATP test to a sensitivity level of 0.1 fmol ATP and results were obtained in 15 seconds to give immediate feedback and corrective action. Surfaces that failed at greater than 10 RLU were re-cleaned and re-tested.

When all surfaces passed with SuperSnap the surfaces were then swabbed with ALLER-Snap protein

Table 2: Cleaning assurance for allergens control using 3 high sensitivity detection methods

Test (LoD)	PRE-CLEANING			POST-CLEANING		
	ATP (RLU) (0.1fmols)	Protein (1 µg)	Sp. Allergen ELISA (16 µg)	ATP (RLU) (0.1fmols)	Protein (1 µg)	Sp. Allergen ELISA (16 µg)
Low risk equipment;	9999	Positive	Positive	2820	Negative	Negative
	Fail			Fail		
	677	Positive	Positive	51	Negative	Negative
	Fail			Fail		
	9999	Positive	Positive	1380	Negative	Negative
	Fail			Fail		
High Care	25	Negative	Negative	17	Negative	Negative
	Pass			Pass		
	2974	Positive	Positive	2	Negative	Negative
	Fail			Pass		
	180	Positive	Positive	11	Negative	Negative
	Fail			Pass		
	2068	Positive	Positive	0	Negative	Negative
	Fail			Pass		
1128	Positive	Positive	3	Negative	Negative	
Fail			Pass			
332	Positive	Positive	12	Negative	Negative	
Fail			Pass			

detection swabs, when they passed, then the specific allergen tests were used. When the protein test gave negative results showing absence at the 1µg level, then the more expensive specific allergen tests were employed. Subsequently, all the specific allergen tests were shown to be negative and the line was released back to production within the 10 days.

The staff felt extremely confident that the specific allergen tests would come back negative following the initial pre-validation using the SuperSnap and ALLER-Snap tests, and may consider pre-releasing before 10 days is up.

Pre-validation screening enabled the site to make significant savings by avoiding repeat testing and further lost production. The Hygiene Manager commented that the combined method approach was very beneficial in releasing the nut production area back into

general production and that: "This process gave me confidence that we would get it right first time with the allergen swabs. This not only saved on cost but more

When all surfaces passed with SuperSnap the surfaces were then swabbed with ALLER-Snap protein detection swabs, when they passed, then the specific allergen tests were used

importantly guaranteed food safety. All our allergen swabs came back clear and the area was released back to general production on plan. I would definitely employ this process again."

The regular use of high sensitivity ATP and high sensitivity protein tests enable high standards of

cleaning to be maintained that can be supplemented with specific allergen tests less frequently and as required.

Summary

Cleaning is one the CCPs for allergen control and a variety of detection methods are available to validate the cleaning processes. Specific allergen tests have their limitations and are expensive whereas other methods have sensitivity but lack specificity.

A combination of three high sensitivity detection methods (ATP, protein and specific allergen tests) provides a more comprehensive sensitivity and rapid result that delivers a timely cost-effective solution.■

Martin Easter
General Manager
Hygiene International Ltd
www.hygiene.net

foodeurope: Keeping you up to date

www.foodmagazine.eu.com

Circulation

We hope you have found this edition of **foodeurope** magazine interesting. It is circulated electronically and as hard copies to shows and conferences, and individual subscribers. We are more than happy to include your colleagues/partners in our electronic distribution list free of charge. Contact john@foodmagazine.eu.com to register your interest.

Subscription and reprints

Alternatively, why not subscribe for hard copies of **foodeurope** magazine? The cost is just €100 per year including postage. We also run a reprint service, enabling you to hand out high quality reprints of your articles to customers and prospects. And it's a lot less expensive than you might think. If you wish to subscribe, or are interested in obtaining reprints, please contact john@foodmagazine.eu.com.

Contributions

We are also interested in receiving good articles for any of our four main sections. If you would like to submit an article, please send an email with your proposed topic, relevant section and a brief summary to

jhoskins@editor.eu.com. Articles should be in English and between 1500–2000 words in length. They should also be original, unpublished work.

Feedback

We are constantly trying to improve **foodeurope** magazine in line with what you the readership want, so if you have any suggestions/ comments/ feedback, please email jhoskins@editor.eu.com. Thank you.

