

For food testing purposes
FOR *IN VITRO* USE ONLY

Resuspension Reagent

Version 1, July 2014

Order No. A 500 14

Reagent for 96 tests

Store the kit at 15 to 25 °C

Table of contents

TABLE OF CONTENTS 2

1. WHAT THIS PRODUCT DOES 3

Contents 3

Storage and Stability 3

2. HOW TO USE THIS PRODUCT 3

3. SUPPLEMENTARY INFORMATION 4

Ordering Information 4

Trademarks 4

Contact and Support 4

4. CHANGE INDEX 4

1. What this Product Does

The Resuspension Reagent is designed for the suspension of microbes in combination with the **foodproof** StarPrep Two Kit. The Reagent eliminates the need for hazardous organic extractions or chaotropic agents.

Contents

Product	Content	Storage
A 500 14	1 container with 35 ml Resuspension Reagent	15 to 25 °C

Storage and Stability

Store the Resuspension Reagent at room temperature (15 – 25 °C) through the expiration date printed on the label.


2. How to Use this Product

The Resuspension Reagent should be used in combination with the **foodproof** StarPrep Two Kit to provide the rapid preparation of bacterial DNA for direct use in PCR.

3. Warranty and Disclaimer of Liability

“Limited Warranty” and “Disclaimer of Liability”: BIOTECON Diagnostics warrants that this product is free from defects in materials and workmanship through the expiration date printed on the label and only if the following are complied with:

- (1) The product is used according to the guidelines and instructions set forth in the product literature;
- (2) BIOTECON Diagnostics does not warrant its product against any and all defects when: the defect is as a result of material or workmanship not provided by BIOTECON Diagnostics; defects caused by misuse or use contrary to the instructions supplied, or if the product is contaminated by improper storage or handling;
- (3) All warranties of merchantability and fitness for a particular purpose, written, oral, expressed or implied, shall extend only for a period of one year from the date of manufacture. There are no other warranties that extend beyond those described on the face of this warranty;
- (4) BIOTECON Diagnostics does not undertake responsibility to any purchaser of its product for any undertaking, representation or warranty made by any dealers or distributors selling its products beyond those herein expressly expressed unless expressed in writing by an officer of BIOTECON Diagnostics;
- (5) BIOTECON Diagnostics does not assume responsibility for incidental or consequential damages, including, but not limited to responsibility for loss of use of this product, removal or replacement labor, loss of time, inconvenience, expense for telephone calls, shipping expenses, loss or damage to property or loss of revenue, personal injuries or wrongful death;
- (6) BIOTECON Diagnostics reserves the right to replace or allow credit for any modules returned under this warranty.



4. Supplementary Information

Ordering Information

BIOTECON Diagnostics is offering a broad range of reagents and services. For a complete overview and for more information, please visit our website at www.bc-diagnostics.com.

Trademarks

foodproof[®] is a trademark of BIOTECON Diagnostics GmbH.

Other brand or product names are trademarks of their respective holders.

Contact and Support

If you have questions or experience problems with this or any other product of BIOTECON Diagnostics, please contact our Technical Support staff (for details see www.bc-diagnostics.com). Our scientists commit themselves to providing rapid and effective help. We also want you to contact us if you have suggestions for enhancing our product performance or using our products in new or specialized ways. Such customer information has repeatedly proven invaluable to us and the worldwide research community.

5. Change Index

First version of the package insert.

BIOTECON Diagnostics GmbH

Hermannswerder 17
14473 Potsdam – Germany
Phone +49 (0) 331 2300-200
Fax +49 (0) 331 2300-299
www.bc-diagnostics.com
bcd@bc-diagnostics.com