


MicroSnap® | *E. coli*

MicroSnap® *E. coli* rapidly detects and enumerates *Escherichia coli* bacteria in food or surface samples. Organisms can be detected in as few as 6 hours, delivering results in the same working day or shift.


Benefits:

- Get results in as little as 6 hours
- Large Enumerable Range
- Real-time data sync and alerts for early risk mitigation
- Data integration with EnSURE® Touch and SureTrend®
- Easy sample collection on liquid, solid and surface tests
- Snap-Valve Technology
- Compact product design takes up minimal lab space
- Reduced cost compared to Petrifilm™ and third-party labs
- Highly sensitive biochemical detection
- Easy-to-read using the highly sensitive EnSURE Touch luminometer
- Optimized workflow with the fewest hands-on steps
- AOAC Research Institute *Performance Tested Method*™ (#071302) validated for raw ground beef, RTE ham*, cooked and raw chicken, raw cod, raw prawns, bagged lettuce*, raw milk, sandwiches, bottled water.
- *Certified for enumeration only
- Store refrigerated (2 to 8 °C)



Sensitivity & Specificity:


- Detects down to <10 CFU in 8 hours or less
- MicroSnap *E. coli* detects *Escherichia coli* bacteria



| Catalog No. | Description | Quantity |
|------------------|--|----------|
| MS2-ECOLI | MicroSnap® <i>E. coli</i> Detection Device | 100 |
| MS1-CEC | MicroSnap® Coliform & <i>E. coli</i> Enrichment Device | 100 |
| MS1-EB-BROTH-9ML | MicroSnap® Enhanced EB Broth | 25 |

Petrifilm® is a trademark of Neogen®

MicroSnap *E. coli* Procedure



See an instructional demo at [HygienaTV](#) or [HygienaTV](#)

Related Products

EnSURE® Touch Monitoring System

EnSURE® Touch is a high-sensitivity luminometer that measures multiple quality tests such as ATP, EB, Total Viable Counts, Coliform, *E. coli*, and more. Plus, send test results data to SureTrend® software to record, trend and analyze results over time. Get real-time notification and take your environmental monitoring program to the next level.



EnSURE® Touch Charging Dock

MicroSnap® Coliform

MicroSnap® Coliform detects and enumerates Coliform bacteria. Results are available in as little as 6 hours. MicroSnap Coliform detects down to 1-10 CFUs.



MicroSnap® Total

MicroSnap® Total delivers total viable organisms test results in 7 hours. MicroSnap Total detects down to <10 CFU of gram-positive and gram-negative aerobic and facultative anaerobic bacteria.



MicroSnap® EB (*Enterobacteriaceae*)

MicroSnap® EB delivers test results in as little as 6 hours. MicroSnap EB detects down to 50 CFUs.



| Catalog No. | Description | Quantity |
|--|--|----------|
| ETOUCH | EnSURE® Touch Monitoring System | 1 |
| www.hygiena.com/suretrend | SureTrend® | Inquire |
| MS1-CEC | MicroSnap® Coliform & <i>E. coli</i> Enrichment Device | 100 |
| MS2-COLIFORM | MicroSnap® Coliform | 100 |
| MS2-ECOLI | MicroSnap® <i>E. coli</i> | 100 |
| MS2-TOTAL | MicroSnap® Total | 100 |
| MS2-EB | MicroSnap® EB | 100 |
| CAL | CalCheck LED Calibration Verification Device | 1 |