

EnSURE™ Touch

Cleaning Verification.
Designed for You.



EnSURE™ Touch

Hygiena understands patient health and maximizing insurance reimbursements are essential to your business. Your success depends on the cleanliness of your facility. Published research shows that even with best efforts put forward, on average, only 34-40% of hospital surfaces are actually cleaned to policy standards. Your constant battle against *C. difficile*, MRSA, and other healthcare-associated infections require a proven and effective cleaning verification program. Implementing a monitoring system improves cleaning thoroughness from 40% to 82%¹. The EnSURE Touch provides you with the data you need to ensure patient health and maximize insurance reimbursements.

Features & Benefits:

EnSURE™ Touch is the next-generation cleaning verification and monitoring system. It collects, analyzes, and reports data to help you prove the effectiveness of your facility cleaning and infection prevention programs. The EnSURE Touch is designed to adapt to your workplace and provides the data you need to optimize infection prevention and lower HAI rates.

- Advanced photodiode technology and unique chemistry for maximum sensitivity and durability
- Intuitive smart navigation makes it easy for anyone to use
- Functions like a smartphone
- Syncs data to SureTrend Cloud
- Responsive 5-inch shatter-proof touch screen
- Built-in tilt sensor for more accurate readings
- Add data fields such as room number to test results
- Supports multiple languages
- Train remote teams with built-in screen sharing capability
- Auto shut off for smart battery management

EnSURE Touch Can Help You

Improve Cleaning

- Measure effectiveness of cleaning tools, products, and procedures
- Identify hard to clean areas

Improve Training

- Provide on-the-spot feedback
- Identify training opportunities

Motivate Teams

- Engage staff in infection prevention
- Promote teamwork among EVS and IP

Track Progress

- Post reports in common areas
- Show performance to management

Reduce HAIs

- Correlate high scores to high HAI risks and low scores to low HAI risks
- Increase reimbursement rate





Innovation at Your Fingertips

Quick test: Run a test in as few as two touches

Locations: Create testing points with customized fields and pass/fail limits on the go without needing desktop software

Plans: Easily group testing locations into plans and randomize test locations to ensure sampling coverage

Users: Manage user access to differentiate personnel from managers. Enhanced administrator control

Results & Reports: Easily view testing results and trend analysis

Search: Quickly find specific testing plans and locations on the fly, reducing time between tests


Sync: Wirelessly synchronize data to SureTrend Cloud

Calibration: Factory-calibrated and preset checks to stay GMP/ISO-compliant

Retest: Swipe through failed test locations and instantly run re-tests to show effectiveness of corrective action

Settings: Customize settings for language, custom fields, test types, and simulate competitor RLU limits

Specifications: Size: 7" x 3" x 1" • Weight: 0.7 lbs • Memory: 2GB internal + cloud storage • Connectivity: USB or Wi-Fi*
Language: Supports multiple languages • Battery: Rechargeable Lithium-Ion battery with USB-C charging *



* Battery life varies based on use and configuration • *USB connectivity coming soon in future software update



The data you need for patient safety and HAI management

SureTrend Cloud is an easy-to-use data analysis software that comes free with every EnSURE Touch. It enables users to monitor testing results and manage user access across multiple facilities, making risk management easier than ever.

Features & Benefits:

- Communicate department performance and achievement to your team and to management
- Fully interactive dashboard allows for quick data viewing, showing you everything at a glance
- Easily manage one or hundreds of instruments from a single SureTrend Cloud account
- Monitor testing across one or multiple facilities
- Full online support and resource center only one click away
- Track re-tests to discover opportunities for retraining or identify equipment that may need replacing
- Quickly create custom reports to trend testing performance
- Automatic updates to SureTrend Cloud ensure you're always running the latest version



suretrend.hygiena.com

EnSURE™ Touch

Enables evidence-based care standards by producing greater cleaning excellence, and...

- Higher HCAHPS top box scores
- Higher CMS star ratings for “How often were the patient’s rooms and bathrooms kept clean?”
- Higher reimbursements from CMS Value-Based Purchasing Program



UltraSnap™ | Surface ATP Test



Catalog No. US2020 | Qty: 100

[▶](#) [▶](#) Instructional Video Online

SuperSnap™ | High-Sensitivity Surface ATP Test



Catalog No. SUS3000 | Qty: 100

[▶](#) [▶](#) Instructional Video Online

AquaSnap™ TOTAL | Water ATP Test



Catalog No. AQ-100X | Qty: 100

[▶](#) [▶](#) Instructional Video Online

AquaSnap™ FREE | Water ATP Test



Catalog No. AQ-100FX | Qty: 100

[▶](#) [▶](#) Instructional Video Online

CalCheck | LED Calibration Verification Device




Catalog No. CAL | Qty: 1

[▶](#) [▶](#) Instructional Video Online

Product Name	Catalog No.	Quantity
EnSURE™ Touch	ETOUCH	1
UltraSnap™ Surface ATP Test	US2020	100
SuperSnap™ High-Sensitivity ATP Surface Test	SUS3000	100
AquaSnap™ Total Water ATP Test	AQ-100X	100
AquaSnap™ Free Water ATP Test	AQ-100FX	100
CalCheck LED Calibration Verification Device	CAL	1

EnSURE Touch is factory calibrated and runs a calibration check on each startup. Recommended calibration control kits available.



Americas:

Hygiena Headquarters

941 Avenida Acaso
Camarillo, CA 93012
Phone: 1-805-388-8007
Fax: 1-805-388-5531
info@hygiena.com

Hygiena Canada

2650 Meadowvale Blvd Unit 14
Mississauga, Ontario L5N 6M5
Phone: 1-833-494-4362 (Toll-free)
or 1-416-686-7962

Hygiena Mexico, S.A. de C.V.

Calle 3 Anegas 409 Bodega 5, Col. Nueva
Industrial Vallejo, Delegación Gustavo A.
Madero, C.P. 07700, CDMX, México.
Phone: +52-5281-4108 or +52-5281-4146

International:

Hygiena International

8, Woodshots Meadow
Watford, Hertfordshire
WD18 8YU, UK
Phone: +44 (0)1923-818821
Fax: +44 (0)1923-818825

Hygiena (Shanghai) Trading Co., Ltd.

Rm.7K, No.518, Shangcheng Rd.
Pudong New District
Shanghai, China 200120
Phone: +86 21-5132-1081, +86 21-5132-1077,
+86 21-5132-1078
Fax: +86 21-5132-1080

Hygiena Diagnóstica España S.L.

P. I. Parque Plata, Calle Cañada Real 31-35,
41900, Camas, Sevilla, Spain
Phone: +34 954-08-1276

help.hygiena.com

 [YouTube.com/HygienaTV](https://www.youtube.com/HygienaTV)

 [Vimeo.com/HygienaTV](https://www.vimeo.com/HygienaTV)

