

Quick and Reliable

Designed specifically for food and beverage analysis labs, GlutenTox[®] Sticks Plus for Reader is a rapid gluten detection kit for food and beverage products. GlutenTox Sticks Plus for Reader contains the G12 anti-gliadin antibody that specifically recognizes the 33-mer peptide, the gluten fraction that triggers celiac disease. Tests are quantified using the **Hygiena Cube Reader**. Each kit contains 12 cassettes.

GlutenTox Sticks Plus for Reader is a practical tool for the control of gluten content in your food products.

GlutenTox Sticks Plus for Reader advantages

- Detection limit 1 ppm (<20 ppm is considered “gluten free”)
- Only one dilution required to quantify
- Quantifies 1-40 ppm of gluten using the Hygiena Cube Reader
- Validated for major food matrices
- Flexible protocols (i.e. polyphenols, hydrolyzed samples, etc.)
- No cross-reactivities
- Results in ≤ 50 minutes for foods and beverages

Industrial applications

- Quantification of gluten presence in foods and beverages
- Gluten-free food quality control
- Trace gluten contamination in food
- Monitoring of gluten-removal
- HACCP, ISO 22000, BRC, and IFS compliance



Product Name	Ordering No.	Reference No.	Quantity
Hygiena Cube	ASY3206	ASY3206	1
GlutenTox Sticks Plus for Reader	KIT3007	KT-6476	12 cassettes