

AlerTox® ELISA Alérgenos

Determinación cuantitativa de alérgenos en **alimentos y bebidas** mediante **kits ELISA Sandwich**

» Alta sensibilidad y especificidad » ELISA Sandwich (96 pocillos)



Kit	LOD	LOQ	Kit	LOD	LOQ
AlerTox ELISA Lupine	0.2 ppm	2 ppm	AlerTox ELISA Pistachio	0.13 ppm	1 ppm
AlerTox ELISA Sesame	0.2 ppm	2 ppm	AlerTox ELISA Egg	0.05 ppm	0.4 ppm
AlerTox ELISA Mustard	1.5 ppm	2 ppm	AlerTox ELISA Lysozyme	2 ppb	25 ppb
AlerTox ELISA Cashew	0.2 ppm	2 ppm	AlerTox ELISA Ovalbumin	5 ppb	25 ppb
AlerTox ELISA Soy (STI)	16 ppb	50 ppb	AlerTox ELISA Casein	0.05 ppm	0.20 ppm
AlerTox ELISA Walnut	0.6 ppm	2 ppm	AlerTox ELISA Milk	0.05 ppm	0.5 ppm
AlerTox ELISA Peanut	0.3 ppm	1 ppm	AlerTox ELISA BLG	1.5 ppb	10 ppb
AlerTox ELISA Hazelnut	0.3 ppm	1 ppm	AlerTox ELISA Crustacean	1 ppb	20 ppb
AlerTox ELISA Almond	0.2 ppm	0.5 ppm	AlerTox ELISA Fish	1.4 ppm	4 ppm
AlerTox ELISA Coconut	0.4 ppm	2 ppm	AlerTox ELISA Macadamia	0.1 ppm	1 ppm

*LOD: Límite de detección
Protocolo corto: 1h30

*LOQ: Límite de cuantificación

Fecha de caducidad: 12 meses desde la fecha de producción