

## Cycler Profile Compatibility for foodproof® Kits

Real-Time PCR Cycler	Profile
<b>Hygiena®</b>	
BAX® System Q7	LP
Dualo 32®	DP
Dualo 32® R <sup>2</sup>	DP
<b>Roche</b>	
LightCycler® 480	LP
LightCycler® 96	LP
<b>Thermo Fisher Scientific</b>	
Applied Biosystems™ 7500 Fast	LP
QuantStudio™ 3 and 5	LP / RP
StepOne™ and StepOnePlus™	LP
<b>Agilent</b>	
AriaMx	LP / DP
Mx3005p	RP
<b>IT-IS</b>	
MyGo Pro	DP
<b>Analytik Jena</b>	
qTower 2.0 / qTower 2.2	RP
qTower <sup>3</sup>	LP / RP
<b>Bio-Rad</b>	
CFX96™	LP
CFX96™ Deep Well	LP / RP

Profile: LP = Low Profile, 0.1 mL, RP = Regular Profile, 0.2 mL, DP = Dualo Profile, 0.1 mL